MECROSS CTO

EN	PTCA Balloon Catheter Instructions for use	Page 2	BG	PTCA балонен катетър Инструкции за употреба	Страница 4
PL	Cewnik balonowy do PTCA Instrukcja użytkowania	Strona 6	FA	بالون کاتتر PTCA دستور العمل استفاده	صفحه 8
DE	PTCA-Ballonkatheter Gebrauchsanweisung	Seite 10	FR	Cathéter à ballonnet pour ACTP Mode d'emploi	Page 12
NL	Gebruiksaanwijzing bij de PTCA-ballonkatheter	Blz. 14	CS	Balónkový katétr pro PTCA Návod k použití	Strana 17
PT	Cateter de balão para PTCA Instruções de utilização	Página 19	ES	Catéter de balón para ACTP Instrucciones de uso	Página 21
EL	Καθετήρας μπαλονιού PTCA Οδηγίες χρήσης	Σελίδα 23	IT	Catetere a palloncino per PTCA Istruzioni per l'uso	Pagina 25

**Legal Manufacturer****Shenzhen Medoo Medical Tech.Co.,Ltd.**

Room A408 Platform Building,
Key Laboratory, Virtual University Park, Yuexing 2 Road,
Yuehai Street, Nanshan District, Shenzhen, Guangdong,
518057, China.

Tel: +86 755-26928966

Fax: +86 755-26920966

Web: www.medoo.hk

**Manufactured in****Jiangsu Medoo Medical Equipment Technology Co.,Ltd.**

No.32 Century Avenue, Dongkan Town, Binhai County, Yancheng,
Jiangsu, 224500, China.

Tel: +86 515-84537776

Fax: +86 515-84537706

Web: www.medoo.hk

**EU Authorized Representative****SUNGO Europe B.V.**

Address : Fascinatio Boulevard 522, Unit 1.7, 2909VA
Capelle aan den IJssel, The Netherlands.

Tel/Fax: +31 (0) 2021 11106

E-mail: ec.rep@sungogroup.com

Web: www.sungogroup.com

2797

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

1.0 Device Description

The McCross CTO (Rx type) device is a coronary dilatation catheter designed for easy guidewire exchange. The catheter working length is 140 cm. Balloon diameters range from 1.0 mm to 4.0mm. The balloon material is made of semi-compliant Pebax material for diameter 1.0 mm to 4.0mm with a rated burst pressure of 16 atmospheres. The proximal shaft of the catheter is composed of a female luer connector bonded to a PTFE coated stainless steel tube with a wire. The proximal shaft joins with a smooth transition to a distal shaft composed of an outer tube and a tri-extrusion inner tube with a balloon laser welded to both tubes at the distal tip. Two radiopaque platinum / iridium marker bands are located within the balloon segment with the exception of balloon diameters less than 2.0 mm which incorporate a centrally positioned single marker band. The inner tube accepts a standard 0.014 inch PTCA guide wire . The guide wire enters the catheter tip and advances coaxially out of the distal Rx port, thereby allowing both coaxial guidance and rapid exchange of catheter with a single standard length guide wire. Two marked sections of 5mm length each located on the proximal shaft indicate catheter position relative to the tip of either a brachial or femoral guiding catheter. The design of this dilatation catheter does not incorporate a lumen for distal dye injections or distal pressure measurements.

2.0 How supplied

Contents	One (1) Balloon Dilatation Catheter One (1) Flushing Needle One (1) Re-wrap tool
Sterile	Sterilized with ethylene oxide gas. Non-pyrogenic
Storage	Store in a dry, dark, normal temperature place.

3.0 Indications

■ The balloon dilatation catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.

4.0 Contraindications

The catheter is contraindicated in:

- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of significant stenosis.

5.0 Warnings

- For single patient, single procedure use only. Do NOT resterilize and/or reuse, as this can potentially result in compromised device performance and increase risk of inappropriate resterilization and cross contamination.
- Do NOT use the catheter if its package has been opened or damaged.
- To reduce the potential for vessel damage in the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA , as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injuries or life-threatening complication.
- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.

- Use the catheter prior to the "Use by" date (Expiration Date) specified on the package.

6.0 Precautions

- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the procedure for which it is to be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty.
- The appropriate anticoagulant and coronary vasodilator should be administered to the patient as needed during the surgical procedure. After the operation, the physician should decide to continue the anticoagulant treatment for a period of time.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn for the patient before removing any additional equipment.
- The used PTCA catheter disposal shall follow individual medical institution/hospital's guidelines.
- Do not reinsert the PTCA catheter into the hoop dispenser after normal use.
- Do not use oil-based contrast agents, organic solvents or alcohol to prevent catheter leakage, damage or lubricant loss.
- The design and construction of the catheter does not provide the user with remote pressure monitoring capabilities.

7.0 Adverse Events

Possible adverse effects include, but are not limited to the following

- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel
- Hemorrhage or hematoma
- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast medium
- Hypo/hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Embolism

8.0 Materials to be used in combination with a balloon catheter include.

- Arterial Sheath
- Femoral or brachial guiding catheter in the appropriate size and configuration
- Hemostatic valve(s)
- Contrast medium diluted 1: 1 with normal saline
- Sterile heparinized normal saline
- 20 cc luer-lock syringe
- Inflation device
- Guide wire diameter not to exceed 0.014"; see product label
- Guide wire introducer
- Guide wire torque device

Note: The above materials are not included with this product.

9.0 Preparation for use

Prior to use, examine all equipment carefully for defects. Examine the dilatation catheter for bends, kinks, or other damage. Do not use any defective equipment. Prepare equipment to be used following manufacturer's instructions or standard procedure. Complete the following steps to prepare the PTCA catheter for use:

1. Remove the protective mandrel from the catheter tip.
2. Slide the protective sheath off the balloon.
3. Flush the guide wire lumen of the PTCA catheter.
4. Attach the syringe with heparinized normal saline to the flushing needle packaged with the catheter; gently insert the needle into the tip of the catheter and flush the guide wire lumen with heparinized normal saline until fluid is seen exiting the guidewire port.
5. Prepare an inflation device with the recommended contrast medium according to the manufacturer's instructions.
6. Evacuate air from the balloon segment using the following procedure.

7. Fill a 20 cc syringe or the inflation device with approximately 4cc of the recommended contrast medium.
8. After attaching the syringe or inflation device to the balloon inflation lumen, orient the dilatation catheter with the distal tip and the balloon pointing in a downward vertical position.
9. Apply negative pressure and aspirate for 15 seconds . Slowly release the pressure to neutral, allowing contrast to fill the shaft of the dilatation catheter.
10. Disconnect the syringe or inflation device from the inflation port of the dilatation catheter.
11. Remove all air from the syringe or inflation device barrel. Reconnect the syringe or inflation device to the inflation port of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
12. Slowly release the device pressure to neutral.
13. Disconnect the 20 cc syringe (if used) and connect the inflation device to inflation port of the dilatation catheter without introducing air into the system.

Caution: All air must be removed from the balloon and displaced with contrast prior to inserting into the body. Otherwise complications may occur.

10.0 Instruction for Use.

1. Insert a guidewire through the hemostatic valve that is on the guiding catheter, following the manufacturer's instructions.
 2. Advance the guide wire carefully into and through the guiding catheter. Withdraw the guide wire introducer, if used.
 3. Attach a torque device to the guide wire, if desired. Under fluoroscopy proceed with accepted PTCA techniques to advance the guide wire to and across the lesion.
 4. Backload the distal tip of the dilatation catheter onto the guide wire ensuring that guide wire exits the catheter at approximately 25cm proximal to the balloon.
 5. Advance the dilatation catheter over the guidewire until it approaches the hemostatic valve.
 6. Open the hemostatic valve, insert the dilatation catheter while maintaining guide wire position and tighten the hemostatic valve. To facilitate insertion , the balloon must be fully deflated to negative pressure.
 7. Tighten the hemostatic valve to create a seal around the dilatation catheter without inhibiting movement of the dilatation catheter. This will allow continuous recording of proximal coronary artery pressure.
- Note: It is important that the hemostatic valve be closed tightly enough to prevent blood leakage around the dilatation catheter shaft, yet not so tight that it restricts the flow of contrast into and out of the balloon or restricts guide wire movement.**
8. Advance the dilatation catheter until the appropriate proximal marker aligns with the hemostatic valve hub. This indicates that the dilatation catheter tip has reached the guiding catheter tip.
 9. Advance the dilatation catheter over the guide wire and into the stenosis. Continue under fluoroscopy and use the radiopaque marker band (s) to position the usable (dilating) section of the balloon within the stenosis.
 10. Continue the procedure using accepted coronary angioplasty technique to dilate the stenosis.

Note: Do not exceed the rated burst pressure printed on the package label. Maintain negative pressure on the balloon between inflations. WARNING: Due to the hydrophilic coating, the balloon may slip away from the target lesion during compression. Care should be taken under high-resolution fluoroscopy so that the balloon does not change position on the target lesion.

11. Withdraw the deflated PTCA catheter and guide wire into the guiding catheter. Using a technique of choice, remove the PTCA catheter, guide wire and guiding catheter from the vasculature. Dispose of the PTCA catheter, guidewire, and guiding catheter appropriately.
12. If the same balloon dilatation catheter is reinserted, the guide wire lumen of the balloon dilatation catheter is flushed using the irrigation needle as described in the Preparation for use section. Before reinsertion, the balloon dilatation catheter should be wiped cleanly with gauze soaked in sterile saline. The balloon can be refolded (after expansion) using the refolding tool provided in the package (attached to the upper right corner of the compliant card). When placing or removing the refolding tool, use the mandrel to support the guidewire lumen and be careful not to damage the balloon. As described in the section on refolding tools.

11.0 Exchange procedure technique

The PTCA catheter has been specifically designed for rapid, single operator balloon exchanges. To perform a dilatation catheter exchange :

1. Loosen the hemostatic valve.
2. Hold the guide wire and hemostatic valve in one hand, while grasping the balloon shaft in the other hand.
3. Maintain guidewire position in the coronary artery by holding the wire stationary, and begin pulling the dilatation catheter out of the guiding catheter while monitoring the wire position under fluoroscopy.
4. Withdraw the deflated dilatation catheter until the guide wire lumen is reached. Carefully pull back the flexible, distal portion of the dilatation catheter out of the rotating hemostatic valve while maintaining the guide wire's position across the lesion.
5. Slide the distal tip of the dilatation catheter out of the hemostatic valve, and tighten valve onto the guide wire to hold it securely in place.
6. Prepare the next dilatation catheter to be used, as previously described in the preparation for use section.
7. Back load another dilatation catheter onto the guidewire as previously described under the instructions for use section, step 4, and continue the procedure accordingly.

12.0 Rewrap Tool

This is an accessory component that allows the balloon to be rewrapped if required.

1. Deflate the balloon by applying negative pressure to the inflation device and maintain under vacuum.
2. Visually inspect the balloon to confirm that it is fully deflated.
3. Removed the Rewrap Tool from Compliance Card.
4. Load the non-flared end of the Rewrap Tool onto the stylet.
5. Carefully load the stylet back through the distal tip of the catheter and past the proximal end of the balloon.
6. While holding the catheter just proximal to the balloon, push the Rewrap device over the balloon in a gentle twisting motion until the entire balloon is covered.
7. Gently remove the rewrap device/stylet assembly.
8. Inspect the balloon for any potential damage. Discard the balloon catheter if there is any visual damage present on the balloon.

13.0 Others

13.1 Packaging requirements

The PTCA balloon dilatation catheter is packaged in sterile packaging. (Oxidized with ethylene oxide), for single use only.

13.2 Reserve & transportation

The product should be stored in a dry, cool place, and the storage environment is well ventilated and free of corrosive gases. The transportation process is strictly protected from heavy pressure, wetness exposure, etc, or transportation requirements as specified in the order contract .

13.3 Sterile

This product is sterilized by ethylene oxide gas and is not pyrophoric. If the package is damaged, do not use the PTCA balloon to dilate the catheter or attempt to sterilize it. Do not use expired products.

13.4 Date of manufacture

See packing.

13.5 Validity period

Sterilization is valid for two years under the specified storage conditions.

14.0 Reference

Physicians should consult recent literature or current medical practice on balloon dilatation, such a published by American College of Cardiology / American Heart Association.

15.0 Disclaimer of Warranty

Descriptions or specifications in Medoo Medical printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Medoo Medical will not be responsible for any direct, incidental or consequential damages resulting from the misuse of the product.

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