


MEDOO 脉动医学







NeuCross RX Ø 1.0-1.5mm distal illustration

ORDERING INFORMATION

Balloon Length (mm)	Balloon Diameter(mm)												
	1.0	1.25	1.5	1.75	2.0	2.25	2.5	2.75	3.0	3.25	3.5	4.0	
8	MC-11-1008	1208	1508	1708	2008	2208	2508	2708	3008	3208	3508	4008	
10	MC-11-1010	1210	1510	1710	2010	2210	2510	2710	3010	3210	3510	4010	
15	MC-11-1015	1215	1515	1715	2015	2215	2515	2715	3015	3215	3515	4015	
20	MC-11-1020	1220	1520	1720	2020	2220	2520	2720	3020	3220	3520	4020	

COMPLIANCE INFORMATION

ø(mm)	1.0	1.25	1.5	1.75	2.0	2.25	2.5	2.75	3.0	3.25	3.5	4.0
pressure(atm)												
2	0.90	1.15	1.40	1.65	1.87	2.15	2.36	2.60	2.82	3.05	3.30	3.78
4	0.95	1.20	1.45	1.70	1.94	2.20	2.45	2.67	2.91	3.15	3.40	3.90
6 (NP)	1.00	1.25	1.50	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	4.00
8	1.05	1.30	1.55	1.80	2.05	2.35	2.57	2.85	3.10	3.35	3.63	4.13
10	1.10	1.35	1.60	1.85	2.11	2.40	2.64	2.95	3.18	3.45	3.74	4.25
12 (RBP)	1.15	1.40	1.65	1.90	2.19	2.50	2.71	3.00	3.27	3.55	3.85	4.35
14	1.20	1.45	1.70	1.95	2.25	2.55	2.80	3.05	3.35	3.65	3.95	4.45



About Medoo

Medoo Medical is a high-tech enterprise who focuses in R&D, production and sales of advanced endovascular interventional medical devices. Founded in 2014, Medoo group is composed of Shenzhen Medoo Medical Tech. Co., Ltd., Shenzhen Medoo Medtech, Co., Ltd., Jiangsu Medoo Medical Equipment Tech. Co., Ltd. and Hongkong Medoo Medical Tech. Co., Ltd.

Medoo Partners have strong technical strength and strong support from Shenzhen Science and Technology Commission, Currently, Medoo has many core projects in the field of vascular intervention, severe heart failure and clinical surgery, including MeCross CTO high-performance semi-compliant pre-dilatation balloon catheter, MeCross NC high-performance non-compliant post-dilatation balloon catheter, MeCross PTA peripheral interventional balloon catheter, NeuCross neuro interventional balloon catheter, MeCross HP high-pressure balloon catheter, TAVI valve balloon catheter, Urinary stone extraction balloon catheter, Esophageal airway dilatation balloon catheter, ultrasonic scalpel, cardiac aid device, PCI accessory, etc.

Integrating the resources of clinical experts in cardiovascular intervention, scientists of Chinese Academy of Science and the R&D specialists on medical consumables of endovascular interventions, Medoo Medical has launched industry-university-research cooperation united with a number of universities and hospitals and established world-class R&D and production centers in Shenzhen and Jiangsu respectively. Medoo Medical has already passed BSI ISO13485 certification (MD: 662153) and GMP certification of China. Its mainstream products have excellent performance, which approved by the European Union CE certification and China NMPA certification, welcomed by the customers all over the world. Medoo Medical not only carries out distribution business in the world with its own brand, but also provides professional OEM/ODM/CDMO production technical services for high-end vascular intervention medical device enterprises at home and abroad.

Culture of Medoo

Pursue excellence in

technology, strive for

perfection and develop

ensure the health & safety create a vigorous Medoo

Strictly abide to laws and Focus on the future life regulations, in order to science and technology,

& mental health, care for human beings, contribute to our society

Treasure the Life

Goals of Medoo

top level products.

- O Develop Medoo Medical to the world's leading national high-end medical device brand
- Benefiting millions of patients with vascular diseases
- O Pursue the harmony and unity of the enterprise, patients and society

Shenzhen Medoo Medical Tech. Co., Ltd.

- Room A-408, Platform Building, Shenzhen Virtual University Park, Yuexing 2 Road, Nanshan district, Shenzhen, China
- +86 755-26928966
- www.medoo.hk

Jiangsu Medoo Medical Equipment Tech. Co., Ltd.

- NO.32 Century Avenue, Donglan Town, Binhai County, Yancheng, Jiangsu, 224500, China
- +86 515-84537776
- zv@medoo.hk

NEUCROSS™RX

RX Intracranial Balloon Dilatation Catheter

NMPA: 国械注准20223031201

Outer S-Coating and Inner Silicone Double Coatings

ensures durable and strong hydrophilic ability to pass through the most complex lesions

> "Seamless" Force Conduction supreme pushability

Small Outer Diameter

Delicate "A-Advance" Tip













RX Intracranial Balloon Dilatation Catheter [NeuCross RX]

Medoo Medical is a world's leading enterprise in the production of both rapid exchange (RX) and over-the-wire (OTW) intracranial balloon dilatation catheters. The NeuCross RX intracranial balloon dilatation catheter (NMPA certification No.: 国械注准20223031201) is suitable for the interventional treatment of intracranial arterial stenosis. Mainly through the femoral artery (or femoral vein, carotid artery, radial artery and other brain related blood vessels) puncture, with the help of digital subtraction angiography system (DSA), under X-ray fluoroscopy, through the puncture equipment and use of special catheter to reach the neck or intracranial blood vessels, diagnosis and treatment of neck or intracranial blood vessels, spinal cord blood vessels lesions.

NeuCross RX has higher security, excellent pushability and passability. In terms of design and materials, Medoo Medical has made a breakthrough. The unique A-Advance tip can quickly pass through the lesion. Super "S-Coating" outer coating and internal silicone double coating make the operation easy. Small folding outer diameter provides extraordinary passability performance. "Seamless" Force Conduction provides supreme pushability. With 48 models (minimum balloon diameter 1.0mm), more specifications, more choices for clinical treatment, highly praised by experts.

Brighten the life. NeuCross RX offers unlimited possibilities for clinical specialists.


(5) TECHNICAL SPECIFICATION

Performance	Specification					
Nominal Pressure	6atm					
Rated Burst Pressure	12atm					
Guide Wire Compatibility	0.014"					
Catheter Compatibility	5F					
	NeuCross RX					

Product Advantages


Higher Security

- RX design, simple operation, faster and more stable.
- Accurate inflation reduces the risk of vascular dissection
- Triple-folding technology, good rewrap performance of 0.030" outer diameter, easier to retreat after expansion, reducing vascular intima damage.
- Suitable with the tortuous vessels, safer to retreat after expansion.





Guidewire lumen oily coating for easier



Extraordinary Crossability

- Ultra-thin balloon with 25% less wall thickness.
- Tip outer diameter 0.016", easily passes through the tortuous vessels and lesions.
- Minimum balloon diameter 1.0 mm, suitable for the most challenging lesions.